

Fast analysis of Dopamine and Serotonin for high time resolution in microdialysis experiments


M. Eysberg¹, L.M. van Heerwaarden², H.-J. Brouwer², N.J. Reinhoud², V. Valentini¹ 1 Antec Scientific, Boston, USA; 2 Antec Scientific, Alphen a/d Rijn, The Netherlands

Improving time resolution for measurements in brain microdialysates

Microdialysis of neurotransmitters in vivo has become an invaluable tool to study neurotransmission in the living brain. Extracellular fluid of the brain is sampled through a the semipermeable membrane with a microdialysis probe.

HPLC analysis requires fractionation of the sample stream, and the size of the fractions will affect time resolution. To accurately measure fast responses, a high time resolution is necessary as show in Fig.1

We developed a robust commercially available on-line solution to improve time resolution, as shown in this poster for DA and 5-HT




On-line sampling of da and 5-ht for analysis in parallel channels


Time resolution is influenced by:

- perfusion rate (typically 1-2 μL/min)
- total analysis time to process a sample
- size and number of serial sample loops

Two sample loops \rightarrow factor 2 better time resolution.



The applied analysis of DA and 5-HT shows reproducible results (<2%RSD) and a detection limit of 100 picomole/L.



Method

The ALEXYS Neurotransmitter Analyzer is equipped with a 10 port valve to simultaneously analyze two different time samples that are collected in two serially installed sample loops.

Conditions for analysis of DA and 5-HT

ALEXYS Neurotransmitter Analyzer with DECADE Elite and 10-port valve 2 mm glassy carbon SenCell, saltbridge reference, AST setting '1' Flow cell 1.5 μL per channel $V_{Injection}$


UHPLC C18 column, 1.0x100 mm, 1.7 µm particles Columns Acetate buffer pH 5.8, ion pairing agent, acetonitrile Mobile phase

175 μL/min at 35 °C; 280 μL/min at 60 °C (pressure 480-490 bar in both cases) Flow rate

Conditions of microdialysis

Perfusion solution 147 mM Na⁺, 1.2 mM Ca²⁺, 3 mM K⁺, 1.2 mM Mg²⁺, 152.4 mM Cl⁻ in water

Temperature influence on analysis time



On-line microdialysis experiment

The tip of the probe had been immersed in perfusion fluid and then transferred for 25 min into standard solution (10 nM DA and 5-HT in perfusion fluid). The response was continuously monitored before, during and after the transfer (Fig. 3).

With a time resolution of 1 data point per 1.8 minutes, fast responses are detected.




Fig. 4 Reponses from sample loop 1 (closed) and 2 (open). The red line indicates the transfer of the

probe from blank to solution containing 10 nM DA and 5-HT.

Conclusion

The ALEXYS UHPLC Neurotransmitter Analyzer is a dedicated system solution to sensitively measure neurotransmitters in small samples.

Sensitivity can now be combined with a time resolution of < 2 min in combination with on-line microdialysis as shown for the analysis of DA and 5-HT.

