Electrochemical Simulation of Phase I Metabolism of Three Novel Cardiovascular Drugs Using UHPLC-MS/MS

<u>Martin Eysberg¹, Małgorzata Szultka-Młyńska², Nico Reinhoud³, Jean-Pierre Chervet³</u> ³Antec Scientific, Zoeterwoude, Netherlands¹Antec Scientific, Boston, MA 02108; ²Department of Environmental Chemistry and Bioanalytics, Nicolaus Copernicus University, 87-100, Poland

Introduction

In this study electrochemistry (EC) coupled with electrospray ionization mass spectrometry (ESI-MS) was used to investigate the metabolic fate of three novel cardiovascular drugs: **Rivaroxaban** (RIV), **Aliskiren** (ALS), and **Prasugrel** (PRS).

Mimicry ("simulation") of the oxidative phase I metabolism was achieved in a simple amperometric thin-layer cell equipped with a boron-doped diamond (Magic Diamond™) working electrode.

Structures of the electrochemically generated metabolites were elucidated from MS/MS experiments. The results obtained by EC-MS were compared with in-vivo studies by analyzing urine samples from patients after administration of the cardiovascular drug.

Instrumental setup & experimental conditions





Figure 1: Schematics instrumental setup for online EC-MS. 1. Syringe pump 2. ROXY Potentiostat

- 3. ReactorCell (zoom)
- 4. Mass Spectrometer

For online EC-MS experiments, the ROXY™ Potentiostat equipped with the ReactorCell™ (Antec Scientific, Boston, USA) and a boron-doped diamond working electrode was used. ALS, RIV, and PRS solutions (100 µM) were prepared in 10mM ammonium formate and adjusted to pH 7.4 with ammonia and acetonitrile (50:50, v/v). The solution was pumped through the electrochemical cell under continuous flow (8 μ L/min) delivered by a syringe pump. The effluent of the cell was introduced into an Agilent 6410 quadrupole electrospray ionization-mass spectrometry (Agilent, Germany), where the positive ionization mode (ESI+) was applied.

Each mass voltammogram was recorded at least three times to ensure reproducibility of the measurements. The mass spectra of the drugs at different applied values of voltage were recorded to observe oxidation processes in real time. Post-processing of the data was performed with the Agilent MassHunter software.

Results 1.1. Electrochemical conversion and optimization






Figure 2: Effect of the working electrode on electrochemical conversion of RIV. MS spectra for: (**A**) RIV without cell (m/z 436) (**B**) MD electrode (**C**) Pt electrode (**D**) GC electrode




Figure 3: Effect of the mobile phase on the signal intensity of target compounds (PRS, RIV, ALS) and their electrochemically generated derivatives/metabolites (A) type of buffer, (B) pH of ammonium formate, and (C) concentration of ammonium formate. **Optimal conditions: 20mM ammonium formate adjusted to pH 7.4**

1.2. Electrochemical oxidation and identification of metabolites

The electrochemical oxidation behavior of the 3 drugs was studied by using the on-line EC/ECI-MS setup. A mass voltammogram was recorded using a continuous voltage ramp from 0 to 3000 mV. The mass voltammogram for RIV is depicted in Figure 4. For ALS and PRS, see reference [1].





Figure 4: (A) Mass voltammogram of RIV. Oxidation of RIV starts at a potential of approximately 1200 mV, evident by the decline of the signal intensity. (B) Main oxidation product ions are formed at voltages > 1600 mV, with m/z 261, 338, 402 and 453. Working electrode: Magic Diamond (BDD); mobile phase: 20mM ammonium formate ((pH=7.4 ± 0.6):acetonitrile (50:50; v/v); flow rate of mobile phase: 8 μL/min. ROXY Potentiostat equipped with ReactorCell (Antec Scientific). MS: Agilent 6410 quadrupole ESI-MS (Agilent,Germany), in positive ionization mode (ESI+).

RIV includes the loss of 16 Da (O), 36 Da (Cl), 150 Da (C4H6NO2, morpholine moiety), and 178 Da (C10H10NO2, morpholine, and benzyl moiety). The main oxidation product, RIV-M5, is related to the hydroxylation reaction, with the formation of a hydroxyl metabolite and probably metabolism occurring at the morpholinone moiety. The formation of the electrochemical metabolite, RIV-M6, presumably involves methylation of the metabolite RIV-M5 and formation of the electrochemical metabolite, RIV-M7, presumably involves simultaneous hydroxylation and methylation of metabolite RIV-M5.









1.3. Analysis of urine samples from patients

Urine samples obtained from patients treated with selected cardiovascular drugs were used for comparison with the proposed electrochemical simulation. After intake of RIV, ALS and PRS, the urine samples were analyzed by UHPLC–MS/MS and the presence of drug and metabolites were confirmed by MS spectra. In case of RIV the drug was detected in the urine samples at a concentration of 5000 ng/mL.

The information gained by the analysis of the urine samples is summarized in Table 1. Four different metabolites of RIV were found in the urine (Fig. 6). Metabolite RIV-M1 was identified as the major metabolite. Metabolites RIV-M1 and RIV-M5 were found in almost equal amounts of 7 to 8% of the dose. Moreover, two other minor metabolites, RIV-M6 and RIV-M7, were observed in the urine samples. But, unchanged rivaroxaban was the major component in the urine. Approximately, 36% of the dose was excreted as unchanged drug in human urine.





Conclusions

In this study, EC–MS was applied as a purely analytical approach to evaluate the oxidation behavior of RIV, ALS, and PRS. We found that:

- On-line EC–MS allows to rapidly generate and identify the oxidation products of the three cardiovascular drugs
- Most metabolites of RIV, ALS and PRS found in-vivo (urine samples) could be generated electrochemically
- Electrochemistry has emerged as an interesting alternative to generate detect and predict the metabolic pathway in early drug development processes
- It can rapidly mimic the main oxidative reactions (Phase 1) that occur in the human body, without biological interactions
- "in-electro" a real alternative to in-vivo and in-vitro drug metabolism

Reference

[1] Małgorzata Szultka-Młyńska; Journal of Chromatography B 1093–1094 (2018) 100–112