


## Introduction

- Polycyclic aromatic hydrocarbons (PAHs) are a group of ubiquitous environmental contaminants
- Responsible for their toxic effects are reactive metabolites: epoxides, radical cations and quinones
- Pyrene was investigated as a PAH model compound (Figure 1)
- Electrochemistry in combination with liquid chromatography and mass spectrometry (EC-LC-MS) was used to study the generation of metabolites and their reactivity towards thiol groups (SH groups)



Pyrene  
Figure 1: Structure of pyrene



Figure 2: Structures of the investigated trapping agents glutathione (GSH) and beta-lactoglobulin A (beta-LGA). Both molecules contain one free thiol group that can act as nucleophile.

## Method




Figure 3: Instrumental on-line setup for the EC-(LC)-MS analysis of pyrene. A) Identification and characterization of oxidation products. B) Studies on reactivity of oxidation products towards SH groups.

- The pyrene containing solution was pumped through an electrochemical thin-layer cell (EC cell) while applying either a potential ramp (recording of mass voltammograms) or a constant potential.
- A) Oxidation products were identified via direct APCI-MS detection (A1) and further characterized by means of HPLC-APCI-MS (A2)
- B) Adduct formation of products was studied by incubation with a trapping agent and subsequent ESI-MS analysis: GSH adducts (B1) and beta-LGA adducts (B2)

## Summary and Conclusion

- EC/MS is a promising tool for simulating metabolic processes of PAHs.
- Quinones were found to be the major oxidation products of pyrene. Distinction between two different isomers was achieved by means of HPLC separation.
- Adduct formation of pyrene quinones with glutathione and beta-lactoglobulin A was observed. This may play a role in the toxicity of PAHs.

## Results and Discussion

### A: Identification and characterization of oxidation products




Figure 4: Mass voltammogram of pyrene. The intensity of the parent compound (m/z 203) decreases, while new signals arise. The major oxidation product with m/z 233 can be attributed to the formation of pyrene quinones (PQ).




Figure 5: Proposed structures of pyrene quinones as known from in vivo studies.<sup>[1]</sup>




Figure 6: HPLC-APCI-MS chromatograms of pyrene before and after oxidation at 1200 mV. The appearance of two different peaks for pyrene quinone (m/z 233) indicates the formation of two isomers.

- Quinones are the major oxidation products of pyrene (Figure 4)
- Two different quinone isomers are formed as can be seen in Figure 6
- Possible structures are shown in Figure 5

### B: Reactivity of pyrene quinones towards thiol groups




Figure 7: Mass voltammogram of pyrene after addition of glutathione (GSH). Besides the signal for pyrene quinone (PQ, m/z 233), additional signals with m/z 538 and m/z 540 arise. Based on exact masses and fragmentation pattern, these signals were identified as glutathione adducts of PQ.




Figure 8: Proposed structures of GSH adducts of pyrene quinone and their formation pathway.




Figure 9: Deconvoluted mass spectra of native beta-LGA and beta-LGA after incubation with pyrene oxidized at 1200 mV. The additional signal with a mass difference of 232 Da compared to the native protein indicates the formation of an adduct of pyrene quinone with beta-LGA.

- Pyrene quinones show high reactivity towards SH groups
- Formation of adducts was observed with glutathione as well as with beta-lactoglobulin (Figures 7 and 9)
- Initial glutathione diol adducts undergo further oxidation to rebuild the quinoid structure (Figure 8) <sup>[2]</sup>

## Literature

- [1] Harper, K., Br J Cancer 1957, 11, 499-507.  
[2] Murty, V. S., Penning, T. M., Chem Biol Interact 1992, 84, 169-188.

## Acknowledgement

We would like to thank the NRW Graduate School of Chemistry (GSC) and the Fonds der Chemischen Industrie (FCI) for financial support and a travel grant for Tina Wigger.