

The finest HPAEC-PAD applications for carbohydrate analysis

Food and Beverage

Mono- and disaccharides
Sugars in meat & fish
Carbohydrates in food according to AOAC
Carbohydrates in instant coffee
Carbohydrates in Dutch candy
Carbohydrates in honey
Oligo- and Polysaccharides in honey
Sugars in beer

Prebiotics Food

Additives
Analysis of Maltodextrin in Syrups
Fructans in infant formula
TGOS in food products
Profiling of FOS

Lactose Free Products

Lactose in dairy & meat
Lactose in lactose-free products

Artificial Sweeteners

Sugar alcohols
Sucralose

Glycoproteins

N-glycans
Monosaccharides
Sialic acids

Sialic acids in glycoproteins

- **ALEXYS™ Carbohydrate Analyzer**
- **SweetSep™ AEX20 anion-exchange column**
- **Acid hydrolysis of sialic acids from glycoproteins**
- **Fast and sensitive analysis of Neu5Ac and Neu5Gc**

Summary

Sialic acids are a class of sugars that plays vital roles in biological functions, including brain development, cell—pathogen interactions, and as disease biomarkers such as cardiovascular diseases and cancer [1,2]. Although commonly found in animal tissues, sialic acids are rarely observed in plants. Over 50 types of sialic acids are known, with N-acetylneuraminic acid (Neu5Ac) and N-glycolyneuraminic acid (Neu5Gc) being the most commonly studied [3]. Given their importance in many biological processes, accurate monitoring Neu5Ac and Neu5Gc levels in glycoproteins is crucial.

In this application note, a HPAEC-PAD method is presented for the compositional analysis of Neu5Ac and Neu5Gc in hydrolyzed glycoprotein samples using the new SweetSep AEX20 anion-exchange column in combination with the ALEXYS carbohydrates analyzer. The method enables fast and sensitive analysis of Neu5Ac and Neu5Gc released from glycoproteins. To demonstrate the applicability of the method, three different glycoproteins were analyzed.

Introduction

Sialic acids are a class of sugars with a nine-carbon backbone. Two of the important sialic acids are N-acetylneuraminic acid (Neu5Ac) and N-glycolylneuraminic acid (Neu5Gc). These sialic acids can be found in many biomolecules such as glycoproteins, glycolipid, or glycosphingolipids [3]. They play crucial roles in stabilizing the protein structure and regulating the ion binding activity and the hydrophilicity [3]. In addition, sialic acids can protect molecules or cells from enzymatic degradation [4, 5].

Due to its important functions, sialic acids are often served as disease biomarkers. For example, elevated amount of sialic acids are associated with metastatic cancer cells [2]. Sialic acids are also part of therapeutic glycoprotein, where they contribute to a better solubility and improved biological activity [5, 6]. It is important to note that humans have a high level of Neu5Gc antibodies, as they typically do not produce Neu5Gc. This can cause immune response when a therapeutic glycoprotein with Neu5Gc is administered [1, 3, 4, 7]. Therefore, accurate analysis of sialic acids in glycoproteins are important, both for diagnostic tool in variety of diseases, and for quality control in therapeutic glycoproteins production.

A novel type of anion-exchange column with SweetSep™ AEX20 stationary phase based on 5 µm polymer particles is used for the compositional analysis of sialic acids released from glycoproteins. Typically, sialic acids are released from glycoproteins using either enzymatic method (neuraminidase) or acid hydrolysis method. In this application note, HPAEC-PAD analysis of sialic acids from glycoproteins with two different acid hydrolysis methods is presented. Several parameters including repeatability, linearity, limits of detection, as well as method's precision and accuracy are evaluated in this application note.

Method

The sialic acid analysis was performed using the ALEXYS™ Carbohydrate Analyzer (Figure 1), consisting of the ET210 eluent tray (for N2 blanketing), a P6.1L quaternary LPG pump, AS6.1L autosampler, CT2.1 column thermostat, and the DECADE Elite electrochemical detector. The SenCell™ with Au working electrode and HyREF (Pd/H₂) reference electrode was selected for sensitive detection of the sialic acids.

Separation

Neu5Ac and Neu5Gc are relatively stronger acids compared to carbohydrates, with a pKa of 2.6 and 2.92, respectively [8]. In



Figure 1. ALEXYS Carbohydrate Analyzer.

an alkaline conditions, Neu5Ac and Neu5Gc are negatively charged, and therefore can be separated using HPAEC. Due to the extreme alkaline conditions, only polymeric anion-exchange columns are suitable for separation of sialic acids. The separation of sialic acids was carried out using SweetSep™ AEX20 columns (4 × 200 mm analytical column and 4 × 50 mm precolumn), which are anion-exchange columns containing highly monodisperse 5 µm ethylvinylbenzene-divinylbenzene copolymer (80% crosslinked) coated with functionalized nanoparticles with dual ion exchange sites (quaternary amine + tertiary amine). In addition, a borate ion trap was installed between the pump and the injector.

The gradient profile described in Table 2 was employed for the separation of the sialic acids. Sialic acids are retained stronger compared to neutral carbohydrates, and therefore required a strong modifier for fast elution [9]. The method started with linear gradient elution at 100 mM NaOH + 70 mM NaOAc from t=0 min to 100 mM NaOH + 300 mM NaOAc at t = 7.5 minute. This mobile phase composition was kept until t = 9 minute to elute any late eluting component and to remove carbonate ions (CO₃²⁻) build up on the column. After the clean-up step the column is equilibrated for 21 minutes to the starting conditions, resulting in a total run time of 30 minutes. The temperature for separation was set at 30°C.

To minimize the introduction of carbonate ions in the mobile phase the eluents were carefully prepared manually using a carbonate-free 50% w/w NaOH solution and electrochemical

Table 1

HPAEC-ECD Conditions

| | |
|--------------------------------|---|
| HPLC system | ALEXYS™ Carbohydrate Analyzer |
| Detector | DECADE™ Elite electrochemical detector |
| Columns | SweetSep™ AEX20, 4 x 200 mm analytical column, 5 μm SweetSep™ AEX20, 4 x 50 mm precolumn, 5 μm Borate ion trap, 4 x 50 mm column, 10 μm (Antec Scientific) |
| Mobile phase (MP) | A: 100 mM NaOH B: 100 mM NaOH + 500 mM NaOAc Eluents blanketed with Nitrogen 5.0 |
| Flow rate | 0.7 mL/min |
| Back pressure | about 200 bar |
| Injection | 10 μL |
| Temperature | 30 °C for separation, 35 °C for detection |
| Flow cell | SenCell with Au WE, stainless steel AE and HyREF (Pd/ H ₂) RE, AST 2 |
| Potential waveform (4-step) | E1, E2, E3, E4: +0.1, -2.0, +0.6, -0.1 V ts, t1, t2, t3, t4: 0.2, 0.4, 0.02, 0.01, 0.07 s |
| I-cell | about 0.2—0.4 μA |
| ADF | 0.1 Hz |
| Range | 2 μA/V |

Table 2

Gradient program

| Time (min) | Mobile phase | %A | %B | Description |
|------------|------------------------------|----|----|---------------------------------------|
| 0 | 100 mM NaOH, 70 mM NaOAc | 86 | 14 | Elution & detection |
| 7.5 | 100 mM NaOH, 300 mM NaOAc | 40 | 60 | |
| 7.5—9 | 100 mM NaOH, 300 mM NaOAc | 40 | 60 | Column clean-up and regeneration |
| 9—30 | 100 mM NaOH, 70 mM NaOAc | 86 | 14 | Equilibration, starting conditions |

grade sodium acetate salt (all commercially available). The diluent was DI water (resistivity >18 MΩ.cm, TOC <5 ppb), which was sparged with Nitrogen 5.0 (purity 99.999%). During analysis the mobile phase headspace is also blanketed with Nitrogen gas (0.2—0.4 bar N₂ overpressure) using the ET210. The inert gas atmosphere will minimize the introduction of CO₂ in the mobile phase and the subsequent formation of CO₃²⁻ ions, ensuring reproducible analysis.

Detection

For the pulsed amperometric detection of the analytes, the Antec SenCell™ electrochemical flow cell is used. This flow cell [10] has a confined wall-jet design and consists of a Au working electrode (WE), HyREF (Pd/H₂) reference electrode (RE) and stainless steel auxiliary electrode (AE). A 4-step potential waveform was applied as described in Table 1. The temperature for detection was set to 35°C. The cell current was typically about 0.2—0.4 μA using these PAD settings under the specified conditions. This particular 4-step waveform with a pulse duration of 500 ms has been claimed to have benefits: (1) a consistent long-term peak area response and (2) minimal electrode wear [11], resulting in less flow cell maintenance and system down time.

Preparation of standards, reagents and samples

Standards: stock standards of Neu5Ac and Neu5Gc were prepared by dissolving 5 mg of the individual compounds in 10 mL H₂O. This corresponds to the concentration of 1.62 mM of Neu5Ac and 1.54 mM of Neu5Gc. Stock standards under these conditions are approximately stable for more than a month in the freezer at -30°C. A combined stock standards in the concentration of 100 μM of Neu5Ac and 10 μM of Neu5Gc were made by mixing 620 μL of Neu5Ac and 65 μM of Neu5Gc and diluting the mix up to 10 mL with H₂O. Working standards were prepared by serial dilution of the combined stock standards with DI water. The working standards were prepared ranging from 0.1 μM—25 μM for Neu5Ac, and 0.01 μM—2.5 μM for Neu5Gc.

Acids for hydrolysis: for releasing sialic acids from the glycoproteins, two different acids were prepared.

Acetic acid, 4M: in a 100 mL volumetric flask, 22.5 mL of glacial acetic acid was diluted with H₂O up to the total volume of 100 mL.

Trifluoroacetic acid (TFA), 0.2%: in a 100 mL volumetric flask, 200 μL of trifluoroacetic acid (99.5%) was diluted with H₂O up to the total volume of 100 mL.

Samples: Three different glycoproteins were subjected to acid hydrolysis to release its sialic acids content:

- Apo-transferrin from bovine
- Fetuin from fetal bovine serum
- α₁-acid glycoprotein from human plasma

The samples were prepared as follows: stock solutions of the individual glycoprotein were prepared by dissolving 3 mg of

glycoprotein in 1 mL water. The glycoprotein stock solutions were frozen in several aliquots to avoid degradation due to freeze-thaw cycle. For hydrolysis, 35 μL of glycoprotein stock solution was mixed with 65 μL H_2O and 100 μL of acid (4M of acetic acid or 0.2% of TFA), resulting in a total volume of 200 μL . The mixture was shortly vortexed and heated in a water bath at 80°C for 2 hours. The samples were dried in a rotary vacuum concentrator at 30°C for 2.5 hours. Note that this hydrolysis method may not be optimal for all samples. The hydrolysis duration as well as delay time between the hydrolysis and evaporation of the acidic medium may affect the sialic acids amount found due to incomplete release or acidic degradation of the free sialic acids [7, 12-14]. The dried samples were reconstituted with 2 mL of H_2O . The samples were then split into two vials, in which one vial was spiked with a known amounts of Neu5Ac and Neu5Gc as described in the Table 3. The samples were ready for injection and analyzed within 24 hours after reconstitution.

Table 3

Added amounts of Neu5Ac and Neu5Gc standards for spike recovery experiment

| Protein | Neu5Ac Added (μM) | Neu5Gc Added (μM) |
|------------------------------|--------------------------------|--------------------------------|
| Apo-transferrin | 1 | 1 |
| Fetuin | 10 | 0.25 |
| α 1-acid glycoprotein | 10 | 1 |

Results

A chromatogram obtained with the 10 μL injection of 10 μM of Neu5Ac and 1 μM of Neu5Gc is depicted in Figure 2. The chromatogram shows baseline separation of Neu5Ac and Neu5Gc. Under this condition Neu5Ac eluted at 4 minute, and Neu5Gc eluted at 6.9 minute. Note, that the use of acetate gradient results in a slight baseline drift.




Figure 2. The chromatogram obtained from a 10 μL injection of 10 μM Neu5Ac and 1 μM of Neu5Gc in DI water.




Figure 3. Calibration curves obtained with 10 μL injections of Neu5Ac (0.1–25 μM) and Neu5Gc (0.01–2.5 μM) standards in DI water.

Table 4

Repeatability of 10 μL injections of Neu5Ac and Neu5Gc standards in DI water (n = 8)

| Compound | RSDs (%) 10 μM Neu5Ac and 1 μM Neu5Gc | | RSDs (%) 1 μM Neu5Ac and 0.1 μM Neu5Gc | |
|----------|---|------|--|------|
| | t_R | Area | t_R | Area |
| Neu5Ac | 0.10 | 0.09 | 0.09 | 0.49 |
| Neu5Gc | 0.06 | 0.21 | 0.06 | 0.78 |

Linearity and repeatability

The linearity was evaluated for the two sialic acid standards in the concentration range of 0.1–25 μM for Neu5Ac and 0.01–2.5 μM for Neu5Gc. Excellent linearity was achieved with correlation coefficients based on the peak area better than 0.999 for both analytes. Calibration curves for Neu5Ac and Neu5Gc are depicted in Figure 3.

The relative standard deviations (RSDs) of the retention time and peak area were determined by repetitive injections of the 10 μM Neu5Ac and 1 μM Neu5Gc standard mix, as well as 1 μM Neu5Ac and 0.1 μM Neu5Gc standard mix in DI water. The retention times were stable as shown in Table 4, with RSD values $\leq 0.11\%$ for both Neu5Ac and Neu5Gc in two different concentrations. The RSDs for peak areas were $< 0.3\%$ for the 10 μM Neu5Ac and 1 μM Neu5Gc standard mix and $< 0.8\%$ for the 1 μM Neu5Ac and 0.1 μM Neu5Gc standard mix. These data demonstrate good reproducibility for the analysis of sialic acids using the presented method.

Limit of detection (LOD) and limit of quantification (LOQ)

The Limit of Detection (LOD) and Limit of Quantification (LOQ) for both sialic acids are shown in Table 5. The LODs were



calculated as the analyte response corresponding to 3× the ASTM noise (average peak-to-peak baseline noise of 10 segments of 0.5 min). The noise was calculated based on a 5-minute section of the baseline from t = 25 minutes to t = 30 minutes. The average response based on 8 replicate injections obtained with the 1 μM Neu5Ac and 0.1 μM Neu5Gc standard mix was used to calculate the LOD and LOQ. The LOQs were calculated in a similar way to LODs, with 10× S/N ratio instead of 3×. The calculated LODs of Neu5Ac (17.0 nM) and Neu5Gc (7.8 nM) demonstrate the excellent sensitivity of the

acetic acid, which is relatively more difficult to remove due to its higher boiling point.

The chromatograms obtained from 10 μL injections of the glycoprotein hydrolysates are shown in Figure 4—9 on the next page. Neu5Ac and Neu5Gc are sufficiently separated from the small, early-eluting peaks in all samples. Both Neu5Ac and Neu5Gc were found in the apo-transferrin and fetuin hydrolysates. As expected, Neu5Gc was not detected in the α₁-acid glycoprotein sample, because the sample is originated from the human plasma. The amount found for Neu5Ac and Neu5Gc from the acetic acid hydrolysates are similar to the corresponding TFA hydrolysates, indicating that a small amount of TFA can be used to release sialic acids from glycoprotein without significant degradation.

The average amount of Neu5Ac and Neu5Gc for all samples are shown in the Table 6. The samples were measured in triplicate and the peak area RSDs for Neu5Ac in the samples range from 0.08% to 1.18%, while for Neu5Gc the peak area RSDs range from 0.34% to 1.36%. It shows that precise measurements of the sialic acids can be achieved with this method. The accuracy of the sialic acid released from glycoprotein cannot be easily evaluated due to the possibility of less optimal hydrolysis. However, comparison with the published value shows that the amount of sialic acids are in the similar range [7]. The method accuracy was assessed based on the sample recovery. The sample recovery can be calculated based on the average amount of the analytes in the sample and spiked sample and compared to the added amount to the spiked sample.

$$\text{Recovery (\%)} = 100\% * \frac{\text{Amount}_{\text{spiked sample}} - \text{Amount}_{\text{sample}}}{\text{Amount}_{\text{standard}}}$$

The sample recoveries are also listed in the Table 6. The sample recovery found for all samples ranged between 90% - 110%.

Table 5

Limit of Detection (LOD) and Limit of Quantitation (LOQ)

| Compound | Limit of detection | | Limit of quantification | |
|----------|--------------------|-----|-------------------------|------|
| | nM | ppb | nM | ppb |
| Neu5Ac | 17.0 | 5.3 | 56.7 | 17.5 |
| Neu5Gc | 7.8 | 2.5 | 25.8 | 8.4 |

method.

Sample analysis

Three glycoprotein samples were hydrolyzed using two different acids. The 2 M acetic acid hydrolysis is a commonly used method to release sialic acids from glycoprotein. It is important to note that in an earlier attempt (data not shown), a small amount of acetic acid left after hydrolysis procedure is detrimental to the released sialic acid, causing appearance of another unknown peak next to the corresponding sialic acid peaks. Hydrolysis using trifluoroacetic acid is a common practice to hydrolyze glycan into its monosaccharides building blocks. However, it is a lesser-known method to release sialic acids from glycoprotein due to its stronger acidic properties [15]. In this application note we described the use of a very low concentration of trifluoroacetic acid as the alternative of using

Table 6

Average sialic acid content in three glycoprotein acid hydrolysates (n = 3)

| Protein | Neu5Ac | | | | Neu5Gc | | | |
|-----------------------------------|--------------------------|--------------|------------------|--------------|--------------------------|--------------|------------------|--------------|
| | Acetic acid hydrolysates | | TFA hydrolysates | | Acetic acid hydrolysates | | TFA hydrolysates | |
| | Amount (μM) | Recovery (%) | Amount (μM) | Recovery (%) | Amount (μM) | Recovery (%) | Amount (μM) | Recovery (%) |
| Apo-transferrin | 0.62 | 110.3 | 0.62 | 108.5 | 0.72 | 102.7 | 0.74 | 101.2 |
| Fetuin | 9.04 | 101.9 | 10.29 | 99.9 | 0.13 | 102.9 | 0.15 | 101.3 |
| α ₁ -acid glycoprotein | 20.53 | 96.2 | 20.58 | 94.3 | n.d. | 103.5 | n.d. | 100.8 |

*n.d. = not detected





Figure 4. The chromatogram of a 10 μ L injection of acetic acid hydrolysates of apo-transferrin sample (black line) and spiked sample (red line).




Figure 7. The chromatogram of a 10 μ L injection of TFA hydrolysates of apo-transferrin sample (black line) and spiked sample (red line).




Figure 5. The chromatogram of a 10 μ L injection of acetic acid hydrolysates of fetuin sample (black line) and spiked sample (red line).




Figure 8. The chromatogram of a 10 μ L injection of TFA hydrolysates of fetuin sample (black line) and spiked sample (red line).




Figure 6. The chromatogram of a 10 μ L injection of acetic acid hydrolysates of α_1 -acid glycoprotein sample (black line) and spiked sample (red line).




Figure 9. The chromatogram of a 10 μ L injection of TFA hydrolysates of α_1 -acid glycoprotein sample (black line) and spiked sample (red line).

References

1. B. Wang, J. Brand-Miller, The role and potential of sialic acid in human nutrition, *European Journal of Clinical Nutrition*, (2003), 57, 1351-1369
2. M.M. Fuster, J.D. Esko, The sweet and sour of cancer: glycans as novel therapeutic targets, *Nature Reviews Cancer*, (2005), 5, 526-542
3. A. L. Lewis, X. Chen, R.L. Schnaar, A. Varki, *Sialic Acids and Other Nonulosonic Acids*, Essential of Glycobiology, Ch.15, 4th edition. Cold Spring Harbor (NY): Cold Spring Harbor Laboratory Press; 2022
4. A. Varki, P. Gagneux, Multifarious roles of sialic acids in immunity, *Annals of The New York Academy of Sciences*, (2012), 1253 (1), 16-36
5. K. Bork, R. Horstkorte, W. Weidemann, Increasing the sialylation of therapeutic glycoproteins: the potential of the sialic acid biosynthetic pathway, *Journal of Pharmaceutical Sciences* (2009), 98 (10), 3499-3508
6. S. Chia, et al., Enhancing pharmacokinetic and pharmacodynamic properties of recombinant therapeutic proteins by manipulation of sialic acid content, *Biomedicine & Pharmacotherapy*, (2023), 163, 114757
7. J. Hu, J. Rohrer, Determination of glycoprotein sialic acid composition using HPAE-PAD in dual eluent generation cartridge mode, ThermoScientific Application note 74121, <https://assets.thermofisher.com/TFS-Assets/CMD/Application-Notes/an-74121-ic-hpae-pad-glycoprotein-sialic-acid-an74121-en.pdf>
8. T. Angata, A. Varki, Chemical diversity in the sialic acids and related α -keto acids: an evolutionary perspective, *Chemical Reviews*, (2002), 102, 439–469
9. Y.C. Lee, Carbohydrate analyses with high-performance anion-exchange chromatography, *Journal of Chromatography A*, (1996), 720, 1-2, 137 – 149
10. Louw, H.J. Brouwer, N. Reinhoud, Electrochemical flowcell, (2016) US patent 9310330
11. R.D. Rocklin, A.P. Clarke, M. Weitzhandler, Improved longterm reproducibility for pulsed amperometric detection of carbohydrates via a new quadruple-potential waveform, *Anal. Chem*, 70 (1998), 1496 – 1501
11. ThermoFisher Scientific, Direct determination of sialic acids in glycoprotein hydrolyzates by HPAE-PAD, AU 180, <https://tools.thermofisher.com/content/sfs/brochures/AU-180-IC-Sialic-Acids-Glycoprotein-Hydrolyzates-AU71730-EN.pdf>
12. W. Zhu, et al., Degradation kinetics and shelf life of N-acetylneuraminic acid at different pH values, *Molecules*, (2020), 25, 5141
13. J.Q. Fan, et al., Comparison of acid hydrolytic conditions for ASN-linked oligosaccharides, *Analytical Biochemistry*, (1994), 219, 375-378
14. F-T.A. Chen, T.S. Dobashi, R.A., Evangelista, Quantitative analysis of sugar constituents of glycoproteins by capillary electrophoresis, *Glycobiology*, (1998), 8 (11), 1045-1052

Conclusion

The ALEXYS™ Carbohydrate Analyzer in combination with the new SweetSep™ AEX20 allows high-resolution separation and sensitive detection of sialic acids released from glycoproteins. Neu5Ac and Neu5Gc are well separated using the presented method and eluted within 7 minutes. The method demonstrates excellent linearity, good repeatability and high detection sensitivity, as evidenced by the low Limit of Detection (LOD) for both sialic acids. Good method accuracy and precision was demonstrated by successful analysis of three different glycoproteins using two different hydrolysis methods. Overall, the presented method enables reliable sialic acid analysis from glycoproteins.

Ordering information

| Detector only | |
|------------------------------------|---|
| 176.0035B | DECADE Elite SCC electrochemical detector |
| 116.4321 | SenCell 2 mm Au HyREF |
| Recommended ALEXYS analyzer | |
| 180.0057W | ALEXYS Carbohydrates Analyzer - gradient (quaternary LPG) |
| 116.4321 | SenCell 2 mm Au HyREF |
| 186.ATC00 | CT2.1 Column Thermostat |
| Column | |
| 260.0020 | SweetSep™ AEX20, 4 x 200 mm column, 5 µm |
| 260.0025 | SweetSep™ AEX20, 4 x 50 mm precolumn, 5 µm |
| 260.0030 | Borate ion trap, 4 x 50 mm column, 10 µm |
| 260.0100* | Pre-column filter PEEK, 0.5 µm |
| Software# | |
| 195.0035 | Clarity CDS single instr. incl LC, AS module |

*) In case samples might contain particulate matter it is advised to use a pre-column filter.

#) Antec ECD drivers are available for Chromeleon CDS, OpenLAB CDS and Empower CDS. The ALEXYS Carbohydrates Analyzer (full system) can also be controlled under Thermo Fisher Scientific Chromeleon™ CDS. Please contact Antec for more details.

For research purpose only. The information shown in this communication is solely to demonstrate the applicability of the ALEXYS system and DECADE Elite detector. The actual performance may be affected by factors beyond Antec's control. Specifications mentioned in this application note are subject to change without further notice.

DECADE Elite, ALEXYS, SenCell, FlexCell and HyREF are trademarks of Antec Scientific. Clarity™ and DataApex™ are trademarks of DataApex Ltd. Chromeleon™ is a trademark of Thermo Fisher Scientific. OpenLAB™ and Chemstation™ are trademarks of Agilent Technologies, Inc. All other trademarks are the property of their respective owners.

Antec Scientific (USA)
 info@AntecScientific.com
 www.AntecScientific.com
 T 888 572 0012

Antec Scientific (worldwide)
 info@AntecScientific.com
 www.AntecScientific.com
 T +31 (172) 268888

