

# Requirements

## for GABA & glutamate analysis

using the ALEXYS<sup>®</sup> Neurotransmitters system (180.0092W) and step gradient option

180.7052u, Edition 5, 2022



## Warning Symbol



The warning sign denotes a warning. It calls attention to a procedure or practice which, if not adhered to, could result in costs, damage or destruction of parts or all of the equipment. Do not proceed beyond a warning sign until the indicated conditions are fully understood and met.

***For research purposes only. The ALEXYS system is not tested by the manufacturer to comply with the In Vitro Diagnostics Directive.***

## Observe safety

Operation of an electrochemical detector can involve the use of hazardous materials including corrosive fluids and flammable liquids. The instrument should only be operated by users with the following expertise:

- Completed degree as chemical laboratory technician or comparable vocational training
- Fundamental knowledge of liquid chromatography
- Knowledge and experience in the safe handling of toxic and corrosive chemicals and knowledge of the application safety measures prescribed for laboratories.
- Participation in an end-user training (daily use of system and chromatography software) performed by the manufacturer or a company authorized by the manufacturer.



**Unskilled, improper, or careless use of the instrument and the related chemicals can create fire hazards, or other hazards which can cause death, serious injury to personnel, or severe damage to equipment and property.**

**Observe all relevant safety practices at all times.**

---

Copyright ©2022, Antec Scientific, The Netherlands. Contents of this publication may not be reproduced in any form or by any means (including electronic storage and retrieval or translation into a foreign language) without prior agreement and written consent from the copyright of the owner. The information contained in this document is subject to change without notice.

DECADE Elite, DECADE Lite, ROXY, ALEXYS, ROXY Exceed, INTRO, FlexCell, ReactorCell,  $\mu$ -PrepCell, ISAAC, HyREF, SenCell and SynthesisCell are trademarks of Antec Scientific. Whatman™ (word and device) and What-man™ (word only) are trademarks of Whatman international Ltd. SOLVENT IFD™ and AQUEOUS IFD™ are trade-marks of Arbor Technologies, Inc. Clarity®, DataApex® are trademarks of DataApex Ltd. Microsoft® and Windows™ are trademarks of Microsoft Corporation. Excel is a registered trademark of the Microsoft Corporation. Chromele-on™ is a trademark of ThermoFisher Scientific. All other trademarks are the property of their respective owners.

The software and the information provided herein is believed to be reliable. Antec Scientific shall not be liable for errors contained herein or for incidental or consequential damages in connection with the furnishing, performance, or use of software or this manual. All use of the software shall be entirely at the user's own risk.

---

# Table of contents

CHAPTER 1 .....	2
Introduction	
CHAPTER 2 .....	3
Chemicals	
Standard stock solutions	3
Microdialysis perfusion solution	3
Calibrator solutions	4
Mobile phase	4
0.1 M borate buffer, pH 10.4	6
1 M sulfite solution	6
OPA reagent	7
System wash solution	8
APPENDIX .....	9

## CHAPTER 1

## Introduction

Thank you for ordering an ALEXYS LC-ECD system. For a successful on-site installation of the GABA-Glu analysis on the ALEXYS system, please arrange the following requirements at your location in advance:

- a computer (see document 195.7000 for the PC requirements)
- general laboratory conditions and facilities, consumables and chemicals for use with the ALEXYS system (see document 180.7070C)
- application specific chemicals and consumables listed in this document



**Arrange these requirements well in advance before the installation to prevent (costly) delays.**

This document lists the application specific chemicals for the analysis of GABA and glutamate using the combination of the following hardware:

ALEXYS® Neurotransmitters system (180.0091U)

ALEXYS GABA/Glu SCC kit (p/n 180.0504)

Step-gradient upgrade kit (p/n 180.0602)

## CHAPTER 2

## Chemicals

For LC-ECD applications, only chemicals of sufficient specific quality should be used to be able to have an optimal system with good performance. The appendix shows detailed descriptions of some of the chemicals that have been used in the Antec R&D laboratory, as an example of what works.



**Have the chemicals and solutions ready at the start of the installation.**

### Standard stock solutions

#### Chemicals

- Standards of GABA and glutamate in high purity grade
- Water (Resistivity >18M $\Omega$ .cm, TOC<10ppb)

#### Preparation

- 10 mL 1 mM GABA in water  
*Store at 4 °C until use (max 1 month)*
- 10 mL 1 mM Glutamate in water  
*Store at 4 °C until use (max 1 month)*

### Microdialysis perfusion solution

In case the sample is in a background of perfusion solution, the standards should also be prepared in this solution.

#### Chemicals

- Chemicals for preparation of microdialysis perfusion fluid
- Water (Resistivity >18M $\Omega$ .cm, TOC<10ppb)

### Preparation

- 100 mL microdialysis perfusion fluid  
Store at 4 °C until use (max 2 month)

Suggestion for perfusion fluid composition:

Component	Concentration (mM)
NaCl	147
KCl	3
MgCl <sub>2</sub> .6H <sub>2</sub> O	1.2
CaCl <sub>2</sub>	1.2

### Calibrator solutions

#### Preparation (on the day of use)

- For a 5-point linear calibration in a relevant concentration range mix and dilute the standard stock solutions in a background of perfusion fluid to the following, 1 mL total volume each:
  - Calibrator 0: Blank (microdialysis perfusion fluid)
  - Calibrator 1: 50 nM GABA; 0.5 uM glutamate
  - Calibrator 2: 100 nM GABA; 1 uM glutamate
  - Calibrator 3: 150 nM GABA; 1.5 uM glutamate
  - Calibrator 4: 200 nM GABA; 2 uM glutamate
  - Calibrator 5: 250 nM GABA; 2.5 uM glutamate

### Mobile phase

#### Mobile phase A (Separation)

#### Mobile phase B (Clean-up)

50 mM phosphoric acid	50 mM phosphoric acid
50 mM citric acid	50 mM citric acid
0,1 mM EDTA	0,1 mM EDTA
pH 3.50	pH 3.50
2% v/v Acetonitrile	50% v/v Acetonitrile

**Note:** these mobile phase compositions are suggested starting conditions and not guaranteed to give best results for all kinds of unknown samples. Condition optimisation may be necessary.

### Chemicals

- Phosphoric acid (commercial solution of 85% w/v in water)
- Citric acid, monohydrate
- Di sodium ethylenediaminetetraacetic acid (EDTA. Na<sub>2</sub>)
- Acetonitrile

- 50% w/w NaOH in water (commercial solution)
- Water (Resistivity >18M $\Omega$ .cm, TOC<10ppb)

### Preparation

#### Base solution

1. Add about 0.9 Litre de-ionized low-TOC water in a wide, clean glass beaker and add a clean stir bar
2. Add 0.0744 g Na<sub>2</sub>EDTA.2 H<sub>2</sub>O and stir till it is completely dissolved
3. Add 6.85 mL 85% w/v phosphoric acid solution
4. Add 21.015 g citric acid and stir till it is completely dissolved.
5. Set the pH of the solution to pH 3.5 using 50% NaOH solution.
6. Add de-ionized low-TOC water to a total volume of 1.0 Litre

*This solution can be stored in the fridge at 4°C for about one week*

#### Mobile phase A

1. Pour 250 mL base solution in a clean glass cylinder
2. Add 10 mL Acetonitrile
3. Add de-ionized low-TOC water to a total volume of 500 mL
4. Transfer the mobile phase to a clean 1L GL45 glass bottle from the OR 110 bottle kit (p/n 184.0202)
5. Degas the mobile phase in a sonic bath for about 15 minutes.

*The mobile phase A should be refreshed at least every 3 days when it is kept at room temperature.*

#### Mobile phase B

1. Take a clean 1L GL45 glass bottle from of the OR 110 bottle kit (p/n 184.0202)
2. Pour 250 mL of the base solution in the 1L bottle
3. Add 250 mL Acetonitrile and stir
4. Degas the mobile phase in a sonic bath for about 15 minutes.

*The mobile phase B should be refreshed at least every week when it is kept at room temperature.*

Check the mobile phases daily for microbial growth. Refresh mobile phase A **at least once every 3 days** or more often if bacterial growth is observed earlier!



**Do not filter the solutions by any means; the 0.2  $\mu$ m inline Whatman filters present in the low pressure solvent lines will take care of filtering.**

0.1 M borate buffer, pH 10.4

#### Chemicals

- Boric acid
- 50% w/w NaOH in water (commercial solution)
- Water (Resistivity >18MΩ.cm, TOC<10ppb)

#### Preparation

1. Add about 90 mL de-ionized low-TOC water in an empty and clean glass beaker and add a clean stir bar
2. Add 0.618 g boric acid and stir till completely dissolved
3. Set the pH of the solution to pH 10.4 using 50% NaOH solution
4. Transfer to a measuring cylinder
5. Add de-ionized low-TOC water to a total volume of 100 mL

*This solution can be stored in the fridge at 4°C for about one month*

1 M sulfite solution

#### Chemicals

- sodium sulfite, anhydrous
- Water (Resistivity >18MΩ.cm, TOC<10ppb)

#### Preparation

1. Weigh 126 mg sodium sulfite and add to an Eppendorf vial
2. Add 1 mL de-ionized low-TOC water
3. Shake vigorously until fully dissolved

*This solution cannot be stored and has to be prepared shortly before preparing the OPA reagent.*

OPA reagent

The derivatization procedure and composition of the OPA reagent was modified from Smith & Sharpe (1994) and Beverly *et al* (2001) resulting in the composition given in Table 1.



Table 1. OPA reagent composition for the analysis of GABA and glutamate.

37 mM OPA 50 mM sodium sulphite 90 mM boric acid, pH 10.4 5% methanol
--

### Chemicals

- Methanol
- Ortho-Phthalaldehyde (OPA).

Note: the amount and composition of impurities in OPA varies from brand to brand. It is strongly recommended to use OPA from Pickering only.

### Preparation

1. Weight 25 mg Ortho-Phthalaldehyde and transfer to a small glass beaker.
2. Add 250  $\mu$ L methanol and shake gently until crystals are dissolved
3. Add 250  $\mu$ L 1M sodium sulfite solution.  
The solution will turn turbid/white.
4. Add 4.5 mL 0.1 M borate buffer, pH 10.4  
The reagent will turn clear again.
5. For the method that uses 9  $\mu$ L of sample, transfer 150  $\mu$ L reagent to a microdialysis/autosampler vial, cap and place in the sample tray compartment.  
For the method that uses only 5  $\mu$ L of sample, the reagent has to be diluted 1:1 with water before placing in the sample tray.

*The OPA reagent is light sensitive. The remaining batch can be covered with aluminum foil and stored for about 2 days in the fridge at 4°C.*

- [1] S. Smith, T. Sharpe (1994) *Journal of Chromatography B*, **652**: 228-233
- [2] L. Beverly, M. G. de Vries, S. D. Bouman, L. M. Arseneau (2001) *Am J Physiol Regulatory Integrative Comp Physiol*, **280**:R563-R569

## System wash solutions

### Chemicals

- Isopropanol
- Acetonitrile
- Water (Resistivity >18MΩ.cm, TOC<10ppb)

### Preparation

- 250 mL water, degassed (autosampler needle wash solution)  
*Cap and store at room temperature until use (max 1 day)*
- 500 mL 20% isopropanol in water, degassed (piston wash)  
*Cap and store at room temperature until use (max 1 month)*
- 250 mL 20% acetonitrile, degassed (column flushing)  
*Cap and store at room temperature until use (max 1 month)*
- 250 mL 50% acetonitrile, degassed (column flushing/storage)  
*Cap and store at room temperature until use (max 1 month)*

## A P P E N D I X

A list of the application specific chemicals with purity and purchase details is shown below as a guideline. The listed brands/purities are not necessarily the best chemicals, but these have been giving good results at the Antec R&D laboratory.

If for any reason alternative chemicals need to be purchased, be aware that chemicals that have a specification of high purity may have been tested for UV-active impurities, which can mean that they may still contain electrochemically active impurities. This is one of the reasons why general 'HPLC grade' water is not suitable for use with EC detection:

- choose chemicals with the same purity or better
- do not choose ultra dry grade chemicals

*Table 2. Brands and purities of chemicals used for application development at Antec.*

Component	Purity	Brand	Order no:	Mw	Kg/L
Ortho-Phosphoric acid, 85% w/v in water	p.a.	Fluka	79620	98.00	D:1.68
Citric acid, monohydrate	p.a.	Acros	124910010	210.14	
Na <sub>2</sub> EDTA. 2 H <sub>2</sub> O	SigmaUltra, 99%	Acros	147855000	372.23	
Acetonitrile	HPLC grade, 99.9%	Acros	268260025	41.05	D:0.781
NaOH, 50% w/v in water	puriss., p.a., for HPLC; 50%	Fluka	71686	40.00	D:1.54
GABA (gamma-aminobutyric acid)	> 99 % (A2129)	Sigma	A2129	103.12	
Glutamate (l-glutamic acid)	> 99% (G1251)	Sigma	49449-100G	147.13	
Ortho-phthalaldehyde (OPA)	97%	Aldrich	P3,940-0	134.13	
OPA (from Pickering)	Chromatographic grade	LC Tech	120	134.13	
sodium sulphite, anhydrous	>98%	Sigma	S,0505	126.0	
boric acid	p.a. (>99.8%)	Merck	1650500	61.83	
Water	TOC <10ppb and deionised, resistivity >18 MOhm-cm (Barnstead Easypure II)				

### Manufacturers/vendors

JT-Baker	<a href="http://www.avantormaterials.com">http://www.avantormaterials.com</a>
Sigma-Aldrich	<a href="http://www.sigmaaldrich.com">http://www.sigmaaldrich.com</a>
Fluka	<a href="http://www.sigmaaldrich.com">http://www.sigmaaldrich.com</a>
Fisher Scientific	<a href="http://www.fishersci.com">http://www.fishersci.com</a>
Barnstead	<a href="http://www.thermoscientific.com">http://www.thermoscientific.com</a>
LC Tech	<a href="http://www.lctech.de">http://www.lctech.de</a>