


Rapid and Sensitive HPAEC-PAD Analysis of Neutral Sugars and Uronic Acids in Biomass

M. Eysberg, C. Marvelous, J. van Schaik, H.J. Brouwer, J.P. Chervet
Antec Scientific, Alphen a/d Rijn, The Netherlands

Introduction



High-carbohydrate biomass source

Polysaccharide extracts

Acid hydrolysis

Sustainable pharmaceuticals
Sustainable energy

Compositional analysis

- ◆ Detailed understanding of biomass carbohydrate composition is an important step in developing its potential application
- ◆ The effectiveness of the extraction process directly influences the biomass's potential applications
- ◆ There is a need for a fast and accurate analysis method to evaluate composition of the extracted carbohydrates.
- ◆ Analysis using High-Performance Anion Exchange Chromatography with Pulsed Amperometric Detection (HPAEC-PAD) in combination with the new SweetSep™ AEX20 anion-exchange column

Stationary Phase

The SweetSep™ AEX20 stationary phase is developed for high-resolution separation of carbohydrates with HPAEC-PAD/MS.

- ◆ Rugged polymeric anion-exchange resin
- ◆ Highly monodisperse latex-coated particles (5 µm)
- ◆ Bifunctional anion exchange sites
- ◆ Fast, high-resolution separation
- ◆ Use of smaller ID column 2.1 mm × 200 mm




Figure 1. SweetSep™ AEX 20 column based on a polymeric stationary phase consisting of monodisperse 5 µm particles coated with latex nano beads, functionalized with quaternary + tertiary amine groups.

Method and instrumentation

HPAEC-ECD conditions

HPAEC system	ALEXYS™ Carbohydrate Analyzer
Columns	SweetSep™ AEX20, 2.1 mm ID × 200 mm analytical column, 5 µm SweetSep™ AEX20, 2.1 mm ID × 50 mm precolumn, 5 µm Borate ion trap, 2.1 mm ID × 50 mm column, 10 µm All columns: Antec Scientific
Mobile phase (MP)	A: 10 mM NaOH B: DI Water (resistivity > 18 MOhm.cm and TOC <5 ppb) C: 200 mM NaOH D: 200 mM NaOAc Eluents blanketed with Nitrogen 5.0
Flow rate	0.18 mL/min
Backpressure	About 190 bar
Injection volume	3 µL
Temperature	27°C for separation, 45°C for detection
Flow cell	SenCell with Au WE, stainless steel AE, and HyREF Pd RE, AST 2
PAD Potential wave-form (4-step)	E1, E2, E3, E4: +0.1, -2.0, +0.6, -0.1 V ts, t1, t2, t3, t4: 0.2, 0.4, 0.02, 0.01, 0.07 s
I-cell	About 0.2—0.5 µA
ADF	0.05 Hz
Range	10 µA/V

Results

Chromatogram of standards




Figure 2. The chromatogram obtained from a 3 µL injection of the 100 µM sugar standard mix in DI water. The black ↑ arrows indicate the start of a baseline elevation. The black ↓ arrows indicate an autozero to remove the baseline current offset. The blue arrow indicates a displacement of hydroxide ions from the column by the eluent containing acetate.

Sample analysis — Acid hydrolysates of wood sample




Figure 3. Chromatograms from 3 µL injections of the acid hydrolysates of wood sample spiked with 5 µM standards (bottom, black line), the acid hydrolysate of wood sample (middle, red line), and the zoom into the baseline of the acid hydrolysate of wood sample (top, blue line).

Sample quantification results (mg/mL) and Limit of Detection (ng/mL)

Compound	Wood	Chitin	LOD	Compound, cont'd	Wood	Chitin	LOD
Mannitol	-	-	3	Fructose	0.003	-	10
Fucose	0.002	-	3	Ribose	-	-	10
Arabinose	0.046	0.026	4	Cellobiose	0.023	-	41
Rhamnose	0.003	-	6	Galacturonic acid	0.003	-	83
Galactose	0.090	-	5	Guluronic acid	-	-	111
Sucrose	-	-	10	Glucuronic acid	-	-	65
Glucose	1.674	-	6	Mannuronic acid	-	-	58
Xylose	0.246	-	5	Iduronic acid	-	-	128
Mannose	0.412	-	7				

Targeted analysis of uronic acids

For applications focusing solely on uronic acids, a fast and targeted analysis using isocratic separation conditions can be applied. This method significantly increases the sample throughput.




Figure 4. Chromatogram obtained from 3 µL injections of the uronic acids standards. Isocratic elution using 34 mM NaOH + 130 mM NaOAc from t = 0 to t = 12 min. LOD values are 43, 62, 31, 37, and 34 ng/mL for galacturonic acid, guluronic acid, glucuronic acid, mannuronic acid, and iduronic acid, respectively.

Conclusion

High-resolution separation of wide range of neutral sugars and uronic acids was achieved within 26 minutes in one chromatographic run. The use of the 2.1 mm ID SweetSep™ AEX20 columns allows for reduced mobile phase consumption and waste, aligning with green analytical chemistry principle. The presented data demonstrate:

- ◆ **Fast, high-resolution separation** of 12 neutral sugars and 5 uronic acids commonly present in biomass
- ◆ **High sensitivity** with limits of detection as low as 3 ng/mL—41 ng/mL (neutral sugars) and 65 ng/mL—128 ng/mL (uronic acids)
- ◆ Extraction efficiency can be optimized based on the sugar composition for further improving biomass utilization
- ◆ The use of ALEXYS Carbohydrate Analyzer in combination with SweetSep AEX20 is the ideal solution for the analysis of biomass hydrolysates

