
## **Complete online electrochemical disulfide bond reduction of antibodies**

## Tomos E. Morgan,\*<sup>a</sup> Craig Jakes,<sup>a,b</sup> Patrice Knightly,<sup>a</sup> Hendrik-Jan Brouwer,<sup>c</sup> Silvia Millán-Martín<sup>a</sup>, Jean-Pierre Chervet,<sup>c</sup> Ken Cook,<sup>d</sup> Sara Carillo,<sup>a</sup> and Jonathan Bones<sup>a,b</sup>

<sup>a</sup>Characterisation and Comparability Laboratory, NIBRT – the National Institute for bioprocessing research and training, Foster Avenue, Mount Merrion, Blackrock, Co. Dublin, A94 X099, Ireland. <sup>b</sup>School of Chemical Engineering and Bioprocessing, University College of Dublin, Belfield, Dublin 4, Ireland. <sup>c</sup>Antec Scientific, Industrieweg 12, 2382NV Zoeterwoude, The Netherlands. <sup>d</sup>Thermo Fisher Scientific, Hemel Hempstead, Herts., UK, HP2 7GE.

Overview—Electrochemical reduction (EC) of antibodies results in cleavage of all inter— and intra-chain disulfide bonds. The online EC reduction included a solvent switching step allowing both EC reduction and LC-MS analysis to be carried out under optimal conditions.



Introduction—Middle-down MS methods require complete reduction of the *inter*- and *intra*-chain disulfide bonds present in antibodies. Reducing the antibody to the subunit level allows easier analytical separation,
HRMS, and tandem MS methods to be used. Electrochemical reduction of disulfide bonds in proteins using a flow-through reactor cell prior to mass spectrometry has been successfully applied as an alternative method to chemical reduction previously.

**Results**—Online EC disulfide bond reduction consistently reduced intact NISTmAb into Lc and Hc subunits. The LC-MS analysis of the antibody could  $\mathbf{n}$ be kept the same as using more traditional chemical reduction techniques as the solvent switching  $\overline{+}$ step allowed separate EC and LC solvents to be optimised.

Intrachain disulfide bond reduction could be easily monitored using exact masses, isotopic, and charge state distributions. Higher charge states, indicative  $\Box$ of increase accessible surface area and relaxed tertiary structure.

An EC potential of 1 V was effective in reducing di-  $\Box$ sulfide bonds, increasing the potential further showed no further increase in disulfide bond reduction but increased the possibility of solvent electrol-**VSIS**.


Interchain disulfide bonds were completely reduced consistently for NISTmAb . However, intrachain disulfide bonds were still intact at room temperature and in aqueous solvents (A, C).


σ


-+


-

Methodology— Integrated online LC-MS EC reduction of Abs. Intact Ab (A) is introduced via a Vanquish Dual Pump F system (B). A discrete EC solvent flows into the ROXY<sup>TM</sup> Exceed potentiostat (C), through the μ-PrepCell SS where EC reduction occurs producing subunits (D). Switching valve setup in 1:6 position (E) traps the subunits onto a short trapping column, before switching positions (F) to backflush the trap with LC-MS compatible solvents onto the LC/MS HRAM Orbitrap<sup>TM</sup> setup (G).









**Conclusion**—Online electrochemical disulfide bond reduction consistently reduced both the *inter*— and *intra*— chain disulfide bonds in a selection of antibodies. The use of a solvent switching valve setup allowed electrochemical reduction to be carried out in the optimal electrochemical conditions followed by LC-MS analysis in the optimal analytical conditions.

References—Nicolardi, S.; et. al., J. Am. Soc. Mass. Spectrom. 2013, 24 (12), 1980-7; Cramer, C. N. et. al., Anal. Chem. 2016, 88 (3), 1585-92. Buter, L. et. al.; J. Chromatogr. A 2017, 1479, 153-160. Switzar, L. et. al; J. Am. Soc. Mass. Spectrom. 2016, 27 (1), 50-8. Zheng, Q. et. al; Int. J. Mass spectrom. 2013, 353, 84-92. Morgan, T., Jakes, C., Brouwer, H-J., Millán-Martín, S., Chervet, J-P., Cook, K., Carillo, S., Bones, J., Inline electrochemical reduction of NISTmAb for middle-up subunit liquid chromatography-mass spectrometry analysis. Analyst, 2021, doi.org/10.1039/D1AN01184G.

