

DISCLAIMER EXTENDED ROXY EC I/O CABLES

The following cable assemblies:

- | | |
|------------------|--|
| 1. P/N 250.0122A | Serial cable, 9M-9F, 5 m straight |
| 2. P/N 188.0501 | Syringe infusion pump serial cable 5 meter |
| 3. P/N 250.0032F | ROXY potentiostat trigger cable 5m |

are non-standard I/O cables for the ROXY EC system and delivered on explicit customer request only.



WARNING: USE AT YOUR OWN RISK

Electromagnetic compatibility: Although the cables passed functional testing (serial communication & I/O functionality) they are not tested/evaluated with respect to EU directive (89/336/EEC) electromagnetic immunity and Emission in combination with the ROXY potentiostat or syringe pump (standard EN61326-1 Electrical equipment for measurement, control and laboratory use - EMC requirements - Part 1: General requirements).

Therefore it is unknown if the extended cables in combination with the intended equipment (potentiostat and syringe pump) are in compliance with the abovementioned directive under all conditions and laboratory environments. In case of EMI interferences please use the serial and I/O cables originally used