

CE Declaration of conformity

We Antec Scientific, Alphen a/d Rijn, The Netherlands, declare that the product:

ET 210 Inert Eluent tray

type 192

to which this declaration relates, is in conformity with the following CE directives:

2015/863/EU Directive on Restriction of the use of certain Hazardous Substances (RoHS)

2012/19/EU Directive on Waste Electrical and Electronic Equipment (WEEE)



Only use the manufacturer-supplied plastic PPCO bottle assemblies (pn 184.0205) in combination with the ET 210 eluent tray. These bottle assemblies are specifically intended for the purpose of inert gas (Helium or Nitrogen) pressurization of the head space above the mobile phase. Before use always check if the bottles are undamaged. **Never apply gas pressure on damaged/scratched bottles.** Manufacturer will not accept any liability for damage, direct or indirect, caused by connecting third-party bottles which do not meet the specifications.

Date: September 21th, 2020

Dr. N.J. Reinhoud (managing Director)