


MEMORY KEY - MANUALS & SOFTWARE

Thanks for your purchase and your confidence in Antec products. We continually strive to minimize our environmental impact from manufacturing, transport and other business activities. In an effort to reduce the use of paper, we are providing all our user manuals and software on a USB Memory key. Please be sure to carefully read the documentation before using the equipment to avoid loss of data, damage to equipment, or personal injury.

The USB memory key contains two main directories 'documents' and 'software' (see below):



For the ROXY product line, all manuals and related documentation are together in the *roxy* directory.

With respect to the available software:

Directory	Software	Which products/procedures apply
clar_8	DataApex Clarity software	DECADE Elite/Lite, ALEXYS LC-EC systems/analyzers
asm	Autosampler Service Manager	AS 110 autosampler
edlg	Dialogue Elite software	DECADE Elite and Lite, ROXY, OQ & PQ
ser	USB-to-serial drivers	HASP license key
chrmln	Chromeleon CDS control driver for DECADE Elite	DECADE Elite/Lite, FW1.11 or higher and ECD driver for Chromeleon driver (pn 176.0105)*
empower	Empower CDS control driver for DECADE Elite	DECADE Elite/Lite, FW1.11 or higher and ECD driver for Empower (pn 176.0110)*
aglnt	Agilent OpenLab CDS control driver for DECADE Elite	DECADE Elite/Lite, FW1.11 or higher and ECD driver for OpenLab (pn 176.0100)*


*) DECADE Elite license for control by this CDS driver has to be purchased separately.

Software and manuals are subjected to changes without notice. For that reason we cannot guarantee that the manuals & software on this USB memory stick are the most recent versions. The most recent versions are available on our website <https://AntecScientific.com>.

DOWNLOAD INSTRUCTIONS

Visit <https://AntecScientific.com> to access the most recent versions of software and manuals:

- The manuals can be downloaded from the section 'Support/Manuals'
- The software can be downloaded from the section 'Products/Software' after registering a (free) account on the Antec web site. You will receive your login credentials by email.



SUPPORT

In case of questions/problems downloading or installing software please contact the Antec support team:

- Per e-mail: support@AntecScientific.com
- Per telephone during office hours:

Antec Scientific, USA

+ 1 (888) 5720012

Antec Scientific, Headquarters, The Netherlands

T +31 172 268888

Antec Scientific 2022. ALEXYS® is a registered trademark of Antec Scientific. Clarity® and DataApex® are trademarks of DataApex Ltd., Microsoft® and Windows™ are trademarks of Microsoft Corporation, OpenLab is a trademark of Agilent technologies, Chromeleon and Xcalibur are trademarks of Thermo Scientific. The software and the information provided herein is believed to be reliable. Antec shall not be liable for errors contained herein or for incidental or consequential damages in connection with the furnishing, performance, or use of software or this manual. All use of the software shall be entirely at the user's own risk.