

WARNING ON USE PRE-SETTABLE TORQUE

The pre-settable torque p/n 250.0064 for assembling of the VT-03 flow cell is a calibrated tool. Recalibration required every 12 months.



Damage to flow cell or any other consequential or indirect damage related to the use of this tool after expiration of calibration date is not covered by warranty.

Copyright ©2010 Antec. All rights reserved. Contents of this publication may not be reproduced in any form or by any means (including electronic storage and retrieval or translation into a foreign language) without prior agreement and written consent from the copyright of the owner. The information contained in this document is subject to change without notice. The information provided herein is believed to be reliable. Antec Leyden shall not be liable for errors contained herein or for incidental or consequential damages in connection with the furnishing, performance, or use of this manual. All use of the hardware or software shall be entirely at the user's own risk.