

# Eluent preparation - A guide to the perfect eluent!

Presenter:	Dr. Detlef Jensen
When:	September 25, 9:30 h
How long?	45 to 60 minutes (max)
Where:	Online Seminar
Registration?	Required
Log in details	It will be provided well in advance.

## Topics of this Webinar

### The participants:

- DI water. A chemical?
  - Why the water quality governs the analytical performance.
- NaOH and CO<sub>2</sub> – an interesting mixture. Why?
  - What can I do about it?
  - Why 50% NaOH?
  - w/w or w/v? A recurring theme!
- Degassing, how and with what?
  - And why is nitrogen suddenly possible?
  - ET210: State-of-the-art in degassing HPAEC-PAD eluents.
- The pushing ion into the solution.
  - What?
  - How and with what, or do it yourself?
  - Filtration
  - and again degassing.
- and even more practical tips.